

Features :~

Touch-key calibration

Automated calibration can be started by pressing keys. (ATX series)
Also, your external calibration weights can be used for span calibration (all models).

Easy setting - Best fit for weighing application

Quickly adjust the desired ratio of stability and response for every application, even during measurement, with one-touch operation.

Expanded Piece Counting function

Unit weights of up to 5 different samples can be easily entered, stored and recalled for use.

Comparator function

Compare samples to target values or pass/fail criteria; results are clearly indicated.

Formulation mode

Convenient for making many measurements of minute samples and seeking the total mass.

Windows Direct Communication Function

Send balance data to Excel® or other Windows® applications without installing any data communication software. By combining standard Auto Print functions with typical spreadsheet functions, even difficult applications can be easily automated. *RS232C interface is needed.



AX



OPERATION MANUAL PLS. CLICK BELOW...

<http://www.shimadzu.com/an/balance/analytical/n9j25k0000015nsz-att/n9j25k0000015w9y.pdf>

Specifications :~

Model	Capacity	Minimum display	Pan size (mm)	Internal calibration	Internal calibration modes
AX200	200g	0.1mg	80dia	✓	Any time with key touch

* Designs and specifications are subject to change for improvement without prior notice.

萬興企業有限公司
BAN HING HOLDING SDN BHD

No. 13-3 Jalan Raja Laut, 50350 Kuala Lumpur, Malaysia
Tel: +603 4048 1000 Fax: +603 4050 8000
WEB: www.banhing.com Email: sales@banhing.com